

Pharmacogenomic Characteristics and IDgenetix-Guided Medication Management for Older Adults with Depression and Anxiety

Feng Cao, PhD, Raymond A. Lorenz, PharmD, BCPP, Andrea E. Hanson, PhD, Robert W. Cook, PhD

Castle Biosciences, Friendswood, TX

Background

- Pharmacogenomic (PGx) testing determines genetic variations in pharmacokinetic and pharmacodynamic genes that impact medication efficacy and tolerability.^{1,2}
- IDgenetix is an advanced 3-in-1 PGx test that incorporates the results of multi-gene testing along with drug-drug interactions and lifestyle factors to optimize medication selection for patients with major depressive disorder (MDD), anxiety, or other neuropsychiatric illnesses.³
- The potential of IDgenetix in older adult populations that are more prone to polypharmacy and adverse drug events has not been fully elucidated.^{4,5}

Methods

- Castle Biosciences laboratory information management system was queried for age, gender, order dates, and ICD-10 diagnosis codes, which included ADHD, anxiety, bipolar disorder, MDD, OCPD, pain, PTSD, and schizophrenia (SCZ).
- Genotypes, phenotypes, and concomitant medications were analyzed for genetic variations, non-genetic interactions, and overall drug recommendations. Drugs were considered actionable if they appeared in the yellow “use with caution and/or increased monitoring” section of the IDgenetix report.

Results

- 737 reports for patients aged 65 and older were analyzed.
- 72% of the patients were female, and the average age was 76 (range 65-103).
- 58% of the patients were on 5 or more medications (average 7).
- Patients averaged 1.6 (range 1-9) ICD-10 diagnosis codes (**Figure 1**).
- All reports in this study had genetic polymorphisms that could impact medication selection in at least 3 of the 15 genes included in the IDgenetix report (**Figures 2-4**).
- Drug-gene interactions accounted for 66% of all drug recommendations, drug-drug interactions accounted for 27%, and lifestyle factors contributed 7% (**Figure 5**).
- The most common antidepressants with actionable recommendations in this population were citalopram and sertraline, while sertraline and escitalopram were the most common concomitant antidepressants (**Table 1**).

Figure 1. ICD-10 Diagnostic Codes




Figure 2. Actionable Genes




Figure 5. Sources of Medication Recommendations




Figure 3. CYP Genetic Polymorphisms





Figure 4. Non-CYP Genetic Polymorphisms



*COMT genotypes: 25% A/A, 45% A/G, 31% G/G; MTHFR A1298C genotypes: 51% A/A, 40% A/C, 9% C/C; MTHFR C677T genotypes: 47% C/C, 41% C/T, 12% T/T; MTHFR enzyme activity: 16% normal, 55% reduced, 29% low

Table 1. The most common antidepressants (left) and the most common concomitant antidepressants (right) with actionable recommendations in the sample population (n=737).

All Antidepressants	Concomitant Antidepressants
citalopram (98%)	sertraline (14%)
sertraline (96%)	escitalopram (9%)
fluoxetine (96%)	duloxetine (8%)
paroxetine (95%)	fluoxetine (4%)
fluvoxamine (95%)	bupropion (4%)

Conclusion

- In patients 65 and older, IDgenetix is being ordered mainly to guide treatment selection for patients with depression and anxiety.
- All patients in this study had multiple genetic polymorphisms potentially impacting medication selection.
- One-third of the IDgenetix-guided medication recommendations were due to drug-drug and lifestyle factor interactions, which is especially important considering the real-world evidence that polypharmacy occurs in almost half of the elderly population.⁶
- Including additional pharmacokinetic and pharmacodynamic genes, beyond CYP2C19 and CYP2D6, along with drug-drug and drug-lifestyle interactions, more than triples the number of clinically relevant medication recommendations.

References

- Brown et al. *Clin Pharmacol Ther.* 2022.
- Swen et al. *The Lancet.* 2023.
- Bradley et al. *J Psychiatr Res.* 2018.
- Saldivar et al. *PGPM.* 2016.
- Sugarman et al. *Drugs Aging.* 2016.
- Wang et al. *Glob Health Res Policy.* 2023.

Disclosures

This study was funded by Castle Biosciences, Inc. FC, RAL, AEH, and RWC are employees and stock/option holders at Castle Biosciences, Inc.