

Continued evaluation of a 31-gene expression profile to predict metastasis in an expanded cohort of 782 cutaneous melanoma patients

Robert W. Cook, PhD, Kyle R. Covington, PhD, Federico A. Monzon, MD
Castle Biosciences, Inc., Friendswood, TX, US

Background

- A 31-gene expression profile (GEP) test that predicts a cutaneous melanoma (CM) tumor's risk of metastasis independent of AJCC stage has been validated as an independent predictor of risk.¹⁻³
- The test provides a binary outcome of low (Class 1) or high (Class 2) risk of metastasis within five years of diagnosis.
- The GEP has been shown to improve staging by identifying over 70% of patients, including clinically and pathologically node-negative patients, who develop distant metastasis or die from their disease.^{1,2}
- We present a multi-center performance study analyzing the GEP in our cumulative cohort of CM patients, and in combination with sentinel lymph node (SLN) status.

Results

Table 1. Cohort demographics

Clinical characteristics	Overall (n = 782)	Class 1 (n=449)	Class 2 (n=333)
Median age (range), years	60 (18-94)	56 (18-94)	64 (19-94)
Development of regional/distant metastasis	224/170	56/40	168/130
Median time to metastasis/follow-up, years	1.3/6.9	2.1/7.4	1.1/5.9
AJCC stage			
I	382 (49%)	326 (85%)	56 (15%)
II	186 (24%)	59 (32%)	127 (68%)
III	207 (26%)	63 (30%)	144 (70%)
IV	7 (1%)	1 (14%)	6 (86%)
Breslow thickness			
Median (range), mm	1.3 (0.1-29.0)	0.7 (0.1-10.0)	2.6 (0.2-29.0)
≤1 mm	318 (40%)	286 (90%)	32 (10%)
>1 mm	458 (59%)	162 (35%)	296 (65%)
unreported	6 (1%)	1 (17%)	5 (83%)
Mitotic index			
<1 mm ²	141 (18%)	111 (79%)	30 (21%)
≥1 mm ²	395 (51%)	176 (45%)	219 (55%)
unreported	246 (31%)	162 (66%)	84 (34%)
Ulceration			
Absent	460 (59%)	332 (72%)	128 (28%)
Present	217 (28%)	44 (20%)	173 (60%)
unreported	105 (13%)	73 (70%)	32 (30%)
Nodal status			
Negative	305 (39%)	159 (52%)	146 (48%)
Positive	201 (26%)	63 (31%)	138 (69%)
unreported	276 (35%)	227 (82%)	49 (18%)

Figure 1. Kaplan-Meier analysis of RFS for GEP predicted outcomes comparing the 782-patient retrospective cohort (A) to three reported prospective studies⁴⁻⁶ (B-D)




Table 2. Cox multivariate regression of survival outcomes in the 782-patient cohort

	RFS			DMFS			MSS		
	HR	95% CI	p value	HR	95% CI	p value	HR	95% CI	p value
Breslow depth	1.2	1.1-1.2	<0.001	1.2	1.1-1.3	<0.001	1.2	1.1-1.4	0.006
Mitotic rate	1.0	1.0-1.0	0.02	1.0	1.0-1.0	0.09	1.0	0.9-1.0	0.6
Ulceration	1.3	0.9-1.9	0.1	1.8	1.2-2.7	0.009	0.9	0.4-1.8	0.8
SLN positive	2.7	1.9-3.7	<0.001	3.1	2.1-4.6	<0.001	3.5	1.8-7.0	<0.001
GEP Class 2	2.7	1.8-4.1	<0.001	2.4	1.5-3.9	<0.001	4.3	1.7-11.0	0.003

Conclusions

- Results achieved in this newly expanded cohort in a multi-center performance study delineate the prognostic accuracy of the 31-gene expression profile test.
- The GEP offers prognostic information that is independent from and adds to conventional staging methods.
- Outcomes in this retrospective cohort are consistent with results from prospective studies of the GEP test.
- GEP result can help guide clinical decision-making regarding risk-appropriate surveillance and follow-up.

References

- Gerami P, Cook RW, Russell MC, et al. *J Am Acad Dermatol* 2015;72(5):780-785.e3.
- Gerami P, Cook RW, Wilkinson J, et al. *Clin Cancer Res* 2015;21(1):175-183
- Zager JS, Messina J, Sondak VK, et al. *J Clin Oncol* 2016;34(suppl); abstr 9581
- Hsueh EC, DeBloom JR, Lee J, et al. *J Hematol Oncol* 2017; in press
- Hsueh EC, Schwartz TL, Lizalek JM, et al. *J Clin Oncol* 2016;34(suppl); abstr 9565
- Greenhaw BN, Brodland DG, Zitelli JA. Proceedings of 48th annual American College of Mohs Surgery Annual Meeting, 2016

Methods

- 261 previously unreported CM primary tumors from 9 centers were analyzed as part of an IRB-approved archival tissue study; data from this cohort were then combined with data from previous validation studies,¹⁻³ exclusive of the training set, for a combined cohort of 782 cases.
- The GEP test was performed in a CAP/CLIA-accredited laboratory using high-throughput RT-PCR assays. Expression data was analyzed using radial basis machine predictive modeling with new cases trained on a previously developed 164-patient cohort.
- Study endpoints included recurrence-free survival (RFS), defined as time to either a regional or distant metastatic event; distant metastasis-free survival (DMFS), defined as time to any metastatic event beyond the regional nodal basin; and melanoma-specific survival (MSS), defined as time from diagnosis to death documented as specifically resulting from melanoma.

Table 3. Kaplan-Meier analysis of survival rates for GEP and SLN predicted outcomes

	RFS		DMFS		MSS	
	5-year rate	# of events	5-year rate	# of events	5-year rate	# of events
Class 1 (n=449)	89%	56 (12%)	93%	40 (9%)	98%	9 (2%)
Class 2 (n=333)	49%	168 (50%)	59%	130 (39%)	81%	51 (15%)
Node negative (n=305)	74%	82 (27%)	83%	56 (18%)	95%	15 (5%)
SLN+ (n=201)	41%	118 (59%)	51%	95 (47%)	77%	39 (19%)

Figure 2. Kaplan-Meier analysis of RFS, DMFS and MSS for SLN and GEP predicted outcomes





Table 4. Net Reclassification Index comparing the addition of GEP to clinical or pathological nodal status

Endpoint	Sensitivity Change	Net Improvement
Recurrence	0.34 (p<0.0001)	0.13 (p<0.0001)
Distant Metastasis	0.32 (p<0.0001)	0.09 (p=0.027)

Figure 3. Schematic representation of risk stratification using AJCC stage with GEP test result to guide patients' clinical management



Acknowledgements

The authors wish to acknowledge Dr. Rene Gonzalez, University of Colorado Cancer Center, Drs. David Lawson, Keith Delman, and Maria Russell, Winship Cancer Institute of Emory University, Dr. Stephen Lyle, University of Massachusetts Medical School, Dr. Gilchrist Jackson, Kelsey-Seybold Clinic, Dr. Anthony Greisinger, Kelsey Research Foundation, Dr. Lee Cranmer, University of Arizona Cancer Center, Dr. T. Christopher Windham, Florida Hospital Memorial Medical Center, Dr. Lewis Kaminester, Dermatology Palm Beach, Dr. Martin Fleming, University of Tennessee Health Science Center, Drs. John Vetto and Sancy Leachman, Oregon Health & Science University, Drs. Laura Ferris and Jonhan Ho, University of Pittsburgh Medical Center, Dr. Alexander Miller, START Center for Cancer Care, Dr. Sarah Estrada, Affiliated Dermatology, Dr. David Pariser, Pariser Dermatology Specialists, Drs. Jane Messina and Jonathan Zager, H. Lee Moffitt Cancer Center, Drs. Jeffrey Wayne and Pedram Gerami, Northwestern University, and Dr. Brian Gastman, Cleveland Clinic for their contribution of samples to the study.

Disclosures

The proprietary GEP test is clinically available through Castle Biosciences as the DecisionDx[®]-Melanoma test (www.SkinMelanoma.com).