

Performance of a diagnostic 35-gene expression profile test (GEP) on difficult-to-diagnose melanocytic lesions

Sarah I. Estrada, M.D.¹, Jeffrey B. Shackelton, M.D.², Howard L. Martin, M.D.³, Brooke H. Russell, Ph.D.⁴, Olga Zolochovska, Ph.D.⁴, Natalie D. Depcik-Smith, M.D.⁵, Nathan J. Cleaver, D.O.⁶, Gregory A. Hosler, M.D., Ph.D.⁷, Matthew S. Goldberg, M.D.^{4,8}
¹Affiliated Dermatology, Scottsdale, AZ. ²Skin Cancer and Dermatology Institute, Reno, NV. ³Sagis, Houston, TX. ⁴Castle Biosciences, Inc., Friendswood, TX. ⁵Aurora Diagnostics GPA Laboratories, Greensboro, NC. ⁶Cleaver Medical Group, Cumming, GA. ⁷ProPath, Dallas, TX. ⁸Icahn School of Medicine at Mount Sinai, NY.

SYNOPSIS

- The accurate diagnosis of melanocytic neoplasms and characterization of their malignant potential remain significant clinical challenges in dermatopathology.¹⁻⁴
- Visual assessment of hematoxylin and eosin (H&E) stained lesions is inherently subjective and relies on expert interpretation and integration of a wide spectrum of architectural and cytologic features that are weighted differently based on the presumed subtype of melanocytic neoplasm and heavily influenced by the pathologists' personal experience and training.⁵
- Difficult-to-diagnose lesions are commonly sent for second opinions to expert dermatopathologists who have more experience with challenging cases; however, the nature of many lesions remains ambiguous with discordant rates of diagnoses ranging from 25-43%.^{1,6}
- The 35-gene expression profile (GEP) test has reported accuracy metrics of 99.1% sensitivity, 96.2% specificity, 96.1% positive predictive value (PPV) and 99.1% negative predictive value (NPV) within the clinically available ≥18-year-old (yo) population (n=474).⁷

OBJECTIVE

To demonstrate use in the intent-to-use difficult-to-diagnose melanocytic lesion population, which includes diagnostically discordant lesions and lesions designated as unknown malignant potential (UMP), cases were analyzed using the 35-GEP test. Accuracy metrics are not provided as a definitive diagnosis for comparison was not achieved.

METHODS

- Samples were acquired from 7 centers. Samples were reviewed by 2 to 3 independent dermatopathologists and used in the study if diagnoses included benign and malignant designations or if diagnoses included more than one unknown malignant potential (UMP) designation.
- Difficult-to-diagnose lesions in the ≥18 yo population which were excluded from analysis from the cohort published in Estrada *et al.* were analyzed with 35-GEP; samples are described in **Table 1**.
- The 35-GEP utilizes dual algorithms based on neural networks to provide a result of benign, intermediate-risk or malignant.⁷ Accuracy metrics cannot be calculated as the true diagnosis was unable to be captured.

RESULTS

Table 1. Demographic information

	Difficult-to-diagnose lesions N=65
Age, median (range)	46 (18-87)
Sex, % male	44.6
Location on body, % (n)	
Abdomen/Chest	6.1 (4)
Acral	3.1 (2)
Back	43.1 (28)
Extremities	32.3 (21)
Head/Neck	13.9 (9)
Other	1.5 (1)

Table 2. 35-GEP results in difficult-to-diagnose lesions

	35-GEP Result		
	Benign, n	Intermediate, n	Malignant, n
Lesions, % (n)	77% (50)	3% (2)	20% (13)

Table 3. Subtype of melanocytic proliferation as determined by submitting dermatopathologist

≥18 yo population	35-GEP Result		
	Benign, n	Intermediate, n	Malignant, n
Blue	4	0	1
Compound	4	0	1
Deep Penetrating	6	0	2
Dysplastic			
Compound	20 ^a	0	1 ^b
Junctional	11 ^c	1	0
Spitz	5	1	8


Dysplastic nevi had different degrees of atypia: a – severe (n=5), moderate (n=8), mild (n=3); b – severe (n=1); c – severe (n=5), moderate (n=3), mild (n=3).

FUNDING & DISCLOSURES


Funding: This study was sponsored by Castle Biosciences, Inc. (CBI), which provided funding to contributing centers for tissue and clinical data retrieval. SIE is a CBI advisor and shareholder. HM, ND-S, NC, GH are CBI advisors. MSG, BHR and OZ are employees and shareholders of CBI.

Figure 1. Representative H&E images


A. Irritated Atypical Spitz




B. Atypical Compound Melanocytic Lesion



C. Compound Dysplastic Melanocytic Nevus



D. Compound Melanocytic Nevus



A. 3 out of 3 UMP concordance, 35-GEP result of intermediate. **B.** 3 out of 3 UMP concordance, 35-GEP result of malignant. **C.** 1 benign and 2 UMP reviews, 35-GEP result of benign. **D.** 3 out of 3 UMP concordance, 35-GEP result of benign.

CONCLUSIONS

- The 35-GEP test is intended to refine diagnoses of melanocytic neoplasms by providing clinicians with an objective ancillary tool with high accuracy.⁷
- In a cohort of 65 difficult-to-diagnose melanocytic lesions, the test provided a result in 100% of lesions, (i.e. no technical failures).
- The 35-GEP test provided an actionable result in 97% of these difficult-to-diagnose lesions, as only 2 cases received an intermediate 35-GEP results.
- Due to the high level of histopathological diagnostic discordance accuracy metrics are not presented.
- Lesion subtypes are presented, however, due to discordance of lesion subtype, only the submitting dermatopathologists' subtype is shown.
- Representative images demonstrate UMPs by histopathological diagnosis that were provided a 35-GEP result.
- A test with these published accuracy metrics in lesions with diagnostic concordance has been shown to alleviate uncertainty in difficult-to-diagnose lesions leading to recommendations for decreased unnecessary procedures while appropriately identifying at-risk patients.⁸

REFERENCES

- Elmore *et al.* *BMJ*. 2017;357:j2813
- Shoo *et al.* *J Am Acad Dermatol*. 2010;62(5):751-756
- Patrawala *et al.* *J Am Acad Dermatol*. 2016;74(1):75-80
- Farmer *et al.* *Hum Pathol*. 1996;27(6):528-531
- Gonzalez *et al.* *J Am Acad Dermatol*. 2017;77(3):543-548
- Piepkorn *et al.* *JAMA Netw Open*. 2019;2(10):e1912597
- Estrada *et al.* *SKIN J Cutan Med*. 2020;4(6):506-522
- Farberg *et al.* *SKIN J Cutan Med*. 2020;4(6):523-533