



Scan or click here for more info

A comprehensive diagnostic offering workflow increases the rate of actionable results of the 23- and 35-gene expression profile tests for use as ancillary diagnostic tools for difficult-to-diagnose melanocytic lesions

Matthew S Goldberg, MD^{1,2}, Jennifer J Siegel, PhD¹, Brooke H Russell, PhD¹, Jason H Rogers, MSc¹, Kyle R Covington, PhD¹, Kristen M Oelschlager, RN¹, Trisha M Poteet¹, Jeffrey K Wilkinson, PhD¹, Michael D Berg, PhD¹, Katherine Falkowski, PhD¹, Sarah J Kurley, PhD¹, and Armineh Kajoian, MD¹

¹Castle Biosciences, Inc., Friendswood, TX ²Icahn School of Medicine at Mount Sinai, New York, NY


Background

- Diagnostic discordance in suspicious cutaneous melanocytic lesions is well documented and particularly prevalent among difficult-to-diagnose cases, for which histopathology may be insufficient for a definitive diagnosis.¹⁻⁴
- The **23-gene expression profile (GEP; myPath Melanoma)** and **35-GEP (DiffDx-Melanoma)** tests are clinically available objective ancillary tools that facilitate diagnosis of melanocytic lesions with ambiguous histopathology. The tests use proprietary algorithms to produce results of likely benign, intermediate, or malignant.⁵⁻⁷
- The 23-GEP has shown accuracy metrics of over 90% sensitivity in multiple clinical studies that included patient outcomes.⁸⁻¹⁰ However, the 23-GEP historically has resulted in ~23% of cases receiving either a technical failure or an intermediate result, which can be perceived as **nonactionable**.^{6,11-13}
- The 35-GEP test can address this shortcoming and showed both an increased sensitivity in the first validation cohort and a decreased nonactionable rate of 8.5%.⁷
- Clinical utility has been demonstrated with benign and malignant GEP test results;^{11,14} therefore, those test results are defined as **actionable**.

Objective

- Today, both the 23- and 35-GEP are offered from a single laboratory as part of a **comprehensive diagnostic offering (CDO)** workflow. Unless preferred otherwise by the ordering clinician, clinical samples are processed first through the 23-GEP test, and if a technical failure or intermediate result is received, processed to the 35-GEP (**Figure 1**).
- Here we report test result metrics from this clinical workflow.

Figure 1. Clinical workflow of the comprehensive diagnostic offering




Methods

- The study included clinical cases submitted to Castle Biosciences for CDO testing with results reported since implementation of the workflow between June 3 and August 31, 2021.
- All cases not receiving a benign or malignant result from the 23-GEP were run on the 35-GEP, except for pediatric cases (<18 years), which were only run on the 23-GEP and excluded from this analysis.
- Technical fail includes samples with insufficient quantity of RNA and/or control or discriminant gene amplification failure based on the requirements for each test.

Results

- The median age for the patient population receiving CDO test results was 49.0 years (range, 18.1-95.5), and 59.2% of patients were female.
- The 23-GEP test gave an actionable result of benign or malignant in 77.8% of cases, which is comparable to past reporting in ambiguous cases for this test^{6,11} (**Figure 2**).
- Nonactionable classifications of the 23-GEP test were 22.2% (12.9% intermediate and 9.4% technical failure). These cases then underwent testing with the 35-GEP test, and an additional 20.9% of originally submitted cases received an actionable result. Only 1.1% of cases received a final intermediate test result (i.e., from both tests); the technical failure rate for the CDO was 0.2% (**Figure 2**).
- This clinical workflow increased the rate of an actionable report from 77.8% to 98.7% when compared with 23-GEP testing alone (**Figure 2**).

Figure 2. Test results from the comprehensive diagnostic offering



Conclusions

- Combining the 23-GEP and 35-GEP tests into a single workflow leverages the strengths of both assays.
- The CDO workflow for histopathologically ambiguous melanocytic lesions has substantially improved reporting of **clinically actionable results** from a historic rate of ~77% for the 23-GEP alone to over 98%.
- Eligible cases with a malignant result from the CDO can also be subsequently run on the 31-GEP prognostic test (Decision Dx-Melanoma), without requiring extra tissue, to predict the likelihood of recurrence and sentinel lymph node biopsy positivity.

References

- Shoo, B. A. et al. *J Am Acad Dermatol* 2010. 62 (5) 751-756.
- Gerami, P. et al. *Am J Surg Pathol* 2010. 34 (6) 816-821.
- Haws, B. et al. *J Cutan Pathol* 2012. 39 (9) 844-849.
- Elmore, J. G. et al. *BMJ* 2017. 357 (1) j2813.
- Clarke, L. E. et al. *J Cutan Pathol* 2015. 42 (4) 244-252.
- Clarke, L. E. et al. *Cancer* 2017. 123 (4) 617-628.
- Estrada, S. et al. *SKIN* 2020. 4 (6) 506-522.
- Ko, J. S. et al. *Cancer Epidem Biomar Prev* 2017. 26 (7) 1107-1113.
- Ko, J. S. et al. *Human Pathology* 2019. 86 213-221.
- Clarke, L. E. et al. *Personalized Medicine* 2020. 17 (5) 361-371.
- Cockerell, C. J. et al. *Medicine* 2016. 95 (40) e4887.
- Minca, E. C. et al. *Mod Pathol* 2016. 29 (8) 832-843.
- Castillo, S. A. et al. *Am J Dermatopathol* 2020. 42 (12) 939-947.
- Farberg, A. et al. *SKIN* 2020. 4 (6) 523-533.

Acknowledgments & Disclosures

- This study was sponsored by Castle Biosciences, Inc.
- All authors are employees and shareholders of Castle Biosciences, Inc.