

Performance of the 35-Gene Expression Profile (GEP) Test for Use as an Adjunctive Melanoma Diagnostic Tool in a New Independent Validation Cohort

Gregory A Hosler, MD, PhD¹, Matthew S Goldberg, MD^{2,3}, Clay J Cockerell, MD⁴, Sarah I Estrada, MD⁵, Howard L Martin, MD⁶, Michael D Berg, PhD², Jennifer J Siegel, PhD², Brooke H Russell, PhD², Sarah J Kurley, PhD², Natalie D Depcik-Smith, MD⁷, Nathan J Cleaver, DO⁸ and Pedram Gerami, MD⁹

¹ProPath, Dallas, TX ²Castle Biosciences, Inc., Friendswood, TX ³Icahn School of Medicine at Mount Sinai, New York, NY ⁴Cockerell Dermatopathology, Dallas, TX ⁵Affiliated Dermatology, Scottsdale, AZ ⁶Sagis, Houston, TX ⁷Aurora Diagnostics GPA Laboratories, Greensboro, NC ⁸Cleaver Medical Group, Cumming, GA ⁹Northwestern University, Chicago, IL


Background

- While histopathological examination is clear for many melanocytic neoplasms, high rates of diagnostic discordance are reported.¹⁻⁴
- Certainty of diagnosis is needed, as the treatment pathways for benign vs. malignant melanocytic neoplasms diverge significantly.
- The **35-gene expression profile (GEP)** test is a clinically available objective ancillary tool that facilitates diagnosis of melanocytic lesions with ambiguous histopathology. The test uses a proprietary algorithm to produce results of likely benign, malignant or intermediate.^{5,6}
- The 35-GEP test showed accuracy metrics of 99.1% sensitivity and 96.4% specificity in the first validation cohort (≥18-years-old, excluding spitzoid lesions).⁵
- The 23-GEP test has shown similar accuracy metrics⁷⁻¹⁰ but has a high test failure rate and a high intermediate result rate.^{8,11,12}
- Both the 35-GEP and the 23-GEP test are now offered through the same laboratory, enabling optimization of clinically actionable results.

Objective

- To describe the performance of the 35-GEP test in a novel cohort independent from those used to develop and initially validate the test⁵ (Figure 1).

Figure 1. Performance Cohort



^aMGF (multiple gene failure) rules: control genes were evaluated independently, and failure of any control gene resulted in sample exclusion. Triplicate gene expression data were aggregated and normalized using control probes.

^bDiagnostically divergent: discordant between benign and malignant designations

^cPediatric: 17 years old or younger

^dSubtype dependent: lesions without a majority subtype or equally split between subtype designation where algorithm inputs change the test result

^eSpitzoid neoplasms: lesions with majority designation of Spitz nevi or spitzoid melanoma.

Methods

- Melanocytic lesions and associated de-identified clinical data from patients ≥18 years of age were obtained from eight laboratories and from lesions submitted for 31-GEP testing (a prognostic melanoma test) in this Institutional Review Board (IRB)-approved study. Lesions were independently reviewed by at least 2 dermatopathologists for diagnostic adjudication and included in the study if they received at least 2 out of 3 diagnostic concordance. Lesion subtype was determined by majority designation and listed as 'subtype discordant' if a majority was not reached.

Results

- Accuracy metrics demonstrate high performance of the 35-GEP test while maintaining a low intermediate result rate (Table 1).
- The 35-GEP demonstrates accurate performance across all subtypes reported, including those with subtype discordance (Table 2).
- The newly modified clinical testing workflow increased the rate of an actionable report from 77.8% to 98.7% when compared with 23-GEP testing alone (Figure 2, Table 3).

Table 1. 35-GEP Performance

	≥18 years old (n=469)	
	35-GEP	95% CI
Sensitivity	96.0%	93.6% - 98.4%
Specificity	92.2%	88.6% - 95.9%
PPV	93.8%	90.9% - 96.7%
NPV	95.0%	91.9% - 98.0%
Intermediate result	2.3%	

Samples that fall in the intermediate zone were excluded from the calculation. PPV, positive predictive value; NPV, negative predictive value; CI, confidence interval.

Table 2. 35-GEP Performance by Subtype

	35-GEP result, n		
	Benign	Intermediate	Malignant
Melanomas			
Acral lentiginous	1	0	25
Desmoplastic	0	0	9
Lentigo maligna	0	0	16
In situ	2	0	18
Nevoid	0	0	2
Nodular	3	0	62
Superficial spreading	1	0	71
Subtype discordant	3	1	39
Nevi			
Blue	33	1	5
Common			
Compound	34	2	0
Intradermal	41	0	1
Junctional	26	3	0
Dysplastic	35 ^a	1 ^b	7 ^c
Subtype discordant	21	3	3

Dysplastic nevi had different degrees of atypia: a - mild (n=1), moderate (n=17) and severe (n=5); b - moderate (n=1); c - mild (n=1), moderate (n=4).

Figure 2. The Comprehensive Diagnostic Offering




Table 3. Actionable Results from Clinical Testing with the Comprehensive Diagnostic Offering

	Actionable (%)	Nonactionable (%)
23-GEP only	77.8%	22.2%
Subsequent 35-GEP	20.8%	1.3%
Overall	98.7%	1.3%

Actionable is the sum of benign and malignant test results. Nonactionable is the sum of intermediate and technical failure test results.

Conclusions

- The 35-GEP test is an accurate, objective ancillary tool that provides information that aids in the diagnosis of melanocytic lesions of uncertain malignant potential.
- This second independent cohort of the 35-GEP test confirms the accuracy metrics compared to the initial validation cohort and demonstrates high accuracy across a range of lesion subtypes.
- The 35-GEP test is now used in conjunction with the 23-GEP test. In clinical cases from June - August, 2021, test results of either benign or malignant were increased from 78% to 99% when the 35-GEP was utilized.
- Eligible cases with a malignant result from the 35-GEP can also be tested on the 31-GEP prognostic test without additional tissue, to predict the likelihood of recurrence and sentinel lymph node biopsy positivity.

References

- Shoo BA et al. *J Am Acad Dermatol*. 2010; 62(5):751-6.
- Gerami P et al. *Am J Surg Pathol*. 2010; 34(6):816-21.
- Haws B et al. *J Cutan Pathol*. 2012; 39(9):844-9.
- Elmore JG et al. *BMJ*. 2017; 357:j2813.
- Estrada S et al. *SKIN*. 2020; 4(6):506-22.
- Farberg A et al. *SKIN*. 2020; 4(6):523-33.
- Clarke LE et al. *J Cutan Pathol*. 2015; 42(4):244-52.
- Clarke LE et al. *Cancer*. 2017; 123(4):617-28.
- Ko JS et al. *Cancer Epidem Biomarker Prev*. 2017; 26(7):1107-13.
- Ko JS et al. *Hum Pathol*. 2019; 86:213-21.
- Minca EC et al. *Mod Pathol*. 2016; 29(8):832-43.
- Castillo SA et al. *Am J Dermatopathol*. 2020; 42(12): 939-47.

Acknowledgments & Disclosures

- This study was sponsored by Castle Biosciences, Inc.
- All authors contributed to and approved the presentation. GAH, CJC, HLM, NDD, NJC, and PG are consultants for Castle Biosciences, Inc. SIE is a consultant and shareholder of Castle Biosciences, Inc. MSG, MDB, JJS, BHR, and SJK are employees and shareholders of Castle Biosciences, Inc.