

# Clinical Utility and a Guide for Using Gene Expression Profile Ancillary Diagnostic Testing for Cutaneous Melanocytic Neoplasms

Joanna Ludzik, MD, PhD<sup>1</sup>, Alexander Witkowski, MD, PhD<sup>1</sup>, Brooke H Russell, PhD<sup>2</sup>, Jason H Rogers, MSc<sup>2</sup>, Matthew S Goldberg, MD<sup>2,3</sup>, Abel D Jarell, MD<sup>4</sup>, and Aaron S Farberg, MD<sup>5</sup>

<sup>1</sup>Department of Dermatology, Oregon Health & Science University, Portland, OR <sup>2</sup>Castle Biosciences, Inc., Friendswood, TX <sup>3</sup>Icahn School of Medicine at Mount Sinai, New York, NY <sup>4</sup>Granite State Dermatology, PC, Portsmouth, NH <sup>5</sup>Baylor Scott & White Health System, Dallas, TX



## Background

- Diagnostic discordance in suspicious cutaneous melanocytic lesions is well documented and particularly prevalent among difficult-to-diagnose cases, for which histopathology may be insufficient for a definitive diagnosis.<sup>1-4</sup>
- The **23-gene expression profile (GEP)** and **35-GEP** tests are clinically available objective ancillary diagnostic tools (collectively, dGEP) that facilitate diagnosis of melanocytic lesions with ambiguous histopathology. The tests use proprietary algorithms to produce results of: **suggestive of benign neoplasm; intermediate (cannot rule out malignancy);** or **suggestive of malignant neoplasm.**<sup>5-7</sup>
- The dGEP tests have demonstrated accuracy metrics of 90.4 - 94.9% sensitivity and 92.5 - 96.2% specificity for 23-GEP, and 94.7 - 99.1% sensitivity and 89.5 - 94.3% specificity for 35-GEP.<sup>5-8</sup>
- Clinical utility has been demonstrated with benign and malignant dGEP test results. For lesions that receive a benign dGEP test result, decreases in treatment intensity have been observed. Specifically, a reduction in the number of excisions performed was noted in two studies (76.7 - 80.5%). A malignant dGEP result prompted 95.2% of dermatologists to increase office visit frequency and a 75% increase in recommendations to excise.<sup>9,10</sup>

## Objective

- Provide a framework for use of dGEP testing in clinical practice
- Demonstrate how GEP results can alter patient management plans

## References


- Shoo, B. A. et al. *J Am Acad Dermatol* 2010. 62 (5) 751-756.
- Gerami, P. et al. *Am J Surg Pathol* 2010. 34 (6) 816-821.
- Haws, B. et al. *J Cutan Pathol* 2012. 39 (9) 844-849.
- Elmore, J. G. et al. *BMJ* 2017. 357 (1)j2813.
- Clarke, L. E. et al. *J Cutan Pathol* 2015. 42 (4) 244-252.
- Clarke, L. E. et al. *Cancer* 2017. 123 (4) 617-628.
- Estrada, S. et al. *SKIN* 2020. 4 (6) 506-522.
- Goldberg, M. et al. *SKIN* 2021.5 s79.
- Cockerell, C. J. et al. *Medicine* 2016. 95 (40) e4887.
- Farberg, A. et al. *SKIN* 2020. 4 (6) 523-533.
- Cockerell C, et al. *Personalized Medicine*. 2017;14(2):123-130.
- Tschen J et al. *Cutis*. 2021; 107(5):264-9.

## Acknowledgments & Disclosures

AW has served as a consultant and investigator for studies supported by Castle Biosciences, Inc. BHR, JHR, and MSG are employees and shareholders of Castle Biosciences, Inc. ASF has served as a consultant for Castle Biosciences, Inc.  
 This study was supported by Castle Biosciences, Inc.

## Results

**Figure 1. A proposed post-biopsy clinical workflow for the use of dGEP for the treating dermatology clinician and the diagnosing dermatopathologist, including situations of diagnostic and clinical use to achieve personalized management and treatment plans.**



**Figure 2. Common provider-based scenarios to consider dGEP testing**

- Dermatopathologist-driven dGEP requests**
  - Atypical melanocytic proliferations
  - Severely dysplastic nevi
  - Atypical blue nevi
  - Clark's or congenital nevi with unusual features
  - Limited tissue availability
- Dermatologist-driven dGEP requests**
  - Significant clinical concern for melanoma e.g., suspicious dermoscopy, confocal, serial imaging
  - Ambiguous pathology report received
  - Personal history of melanoma
  - Cosmetically sensitive areas

**Figure 3. Overall Summary of dGEP Diagnostic and Treatment Utility**

Diagnostic Utility	Treatment Utility
Is it melanoma or not?	Should an excision be performed?
<ul style="list-style-type: none"> <li>43% reduction in indeterminate diagnoses<sup>11</sup></li> <li>51% increase in diagnostic confidence<sup>10</sup></li> <li>High accuracy<sup>5,8</sup></li> </ul>	<p><b>Benign dGEP Result</b></p> <ul style="list-style-type: none"> <li>76.7%-80.5% reduction in excisions<sup>9,10</sup></li> <li>74.1% of dermatologists reduce office visits<sup>10</sup></li> <li>Demonstrated safe to forego re-excision<sup>12</sup></li> </ul> <p><b>Malignant dGEP Result</b></p> <ul style="list-style-type: none"> <li>75% increase in excisions<sup>9,10</sup></li> <li>95.2% of dermatologists increase office visits<sup>10</sup></li> </ul>

## Conclusions

- GEP testing and results can be incorporated into diagnostic clinical use at various points during the patient journey.
- Optimal dGEP use takes place during the diagnostic process and involves the dermatopathologist to alleviate ambiguous diagnoses.
- Clinical concerns that raise uncertainty may necessitate utilization of dGEP to achieve clinicopathological correlation.
- Treatment decisions have been shown to be dGEP result-dependent and can alter patient management plans such as excision and office visit frequency.