

# The 40-gene expression profile (40-GEP) test allows for an improved prognostication of the likelihood of metastasis in patients with T1 cutaneous squamous cell carcinoma (cSCC) with high-risk factors

Ally-Khan Somani, MD, PhD<sup>1</sup>; Javier Cañueto MD, PhD<sup>2</sup>; Vincent Smith, DNP, ANP-C<sup>3</sup>; Andrew Berg, NP-C<sup>4</sup>; Alison L. Fitzgerald, PhD<sup>5</sup>; Jennifer J. Siegel, PhD<sup>5</sup>; Anesh Prasai, PhD<sup>5</sup>; Sherrif F. Ibrahim, MD, PhD<sup>6</sup>; Sarah T. Arron, MD, PhD<sup>7</sup>

<sup>1</sup> Department of Dermatology, Indiana University School of Medicine, Indianapolis, Indiana, USA; <sup>2</sup> Department of Dermatology, Universitario de Salamanca, Salamanca, Spain; <sup>3</sup> Pinnacle Dermatology, Little Rock, Arkansas, USA; <sup>4</sup> Dermatology Specialists of Kansas City, Kansas City, MO, USA; <sup>5</sup> Research and Development, Castle Biosciences, Inc., Friendswood, Texas, USA; <sup>6</sup> Rochester Dermatologic Surgery, Victor, NY, USA; <sup>7</sup> Peninsula Dermatology, Burlingame, California, USA

## Synopsis

- › The metastatic rate for cSCC is low, however, the overall incidence is high (~1-2.5 million cases/year), and deaths from this disease now surpass those from melanoma.<sup>1,2</sup>
- › Guidelines recommend that patient management decisions should be determined by the individual patient's probability of disease progression. However existing prognostic tools are variable and have demonstrated low accuracy, thus a significant number of risky patients are missed.<sup>3-5</sup> (Figure 1)
- › The 40-gene expression profile (40-GEP) test stratifies primary cSCC patients having one or more clinicopathologic risk factors into three biological risk groups based on risk for regional, nodal, or distant metastasis 3-years post diagnosis (low = Class 1; moderate = Class 2A; high = Class 2B).<sup>6</sup> (Figure 2)
- › Studies have demonstrated that when the 40-GEP test is incorporated with staging, it enhances risk assessment and offers a more accurate and personalized stratification of metastatic risk.<sup>6-8</sup>


**Figure 1. Variability in risk factor assessment for cSCC impacts staging and therefore treatment decisions**

Clinicopathologic risk factor	40-GEP testing criteria	Factors used for risk assessment*		
		NCCN (v1.2023)	AJCC8	BWH
Tumor size ≥2 cm	✓	✓	✓	✓
Invasion beyond subcutaneous fat or >6mm**	✓	✓	✓	✓
Perineural invasion#	✓	✓	✓	✓
Poorly differentiated	✓	✓	✓	✓
Recurrent†	✓	✓	✓	✓
Immunosuppression	✓	✓	✓	✓
Site of prior RT or chronic inflammation	✓	✓	✓	✓
Located on head, neck, anogenital, hands, and feet, any size	✓	✓	✓	✓
Borders poorly defined	✓	✓	✓	✓
Rapidly growing tumor	✓	✓	✓	✓
Neurological symptoms	✓	✓	✓	✓
Lymphatic or vascular involvement	✓	✓	✓	✓
Desmoplastic SCC	✓	✓	✓	✓
Specific high-risk subtypes##	✓	✓	✓	✓

**Blue shaded area connotes T1 Stage Tumors**

\*NCCN: National Comprehensive Cancer Network guidelines (high or very-high risk factors); AJCC8: American Joint Committee on Cancer Eighth Edition; BWH: Brigham and Women's Hospital \*\*AJCC8 and NCCN: >6mm and bone erosion/invasion included; BWH: bone invasion automatically upstages to highest risk stage- T3; #AJCC8 and NCCN: ≥0.1mm nerve or deeper than dermis; BWH: ≥0.1mm nerve required; ##Acantholytic, adenosquamous, or metaplastic subtypes (40-GEP: others will be considered on a case-by-case basis); †40-GEP is not validated for local recurrence

**Figure 2. Clinical use of the 40-GEP test**



## Methods


- › Retrospective study consisted of cSCC patients with confirmed clinicopathologic information, which met requirements for 40-GEP clinical testing and follow up of at least 3 years (n=954). 446 cases were staged as BWH T1 (see Figure 1) and 501 cases were staged as AJCC8 T1 (Figure 1).
- › Clinicopathologic risk factors were comprehensively assessed by performing the following on the T1 subsets: review of original biopsy and surgical reports, independent review by a board-certified dermatopathologist, and Kaplan-Meier analysis for metastasis-free survival (MFS) with log rank test and Cox regression analyses.
- › For real world clinical usage data, summary metrics on T1 (n=2,924) samples were identified from all n=7317 patients who received test reports during the first two years of clinical ordering (August 31, 2020-August 31, 2022).

## Clinical Issue and Objective

Staging systems for cSCC have shown limited accuracy in predicting the prognosis of the disease to metastasize.<sup>3,4</sup> Approximately 30% of BWH or AJCC8 T1 cSCCs (deemed at a low risk of metastasis) will go on to metastasize, potentially due to risk factors not captured in clinicopathologic risk assessment, but within the intended use of 40-GEP.<sup>9,10</sup> The lack of accurate risk assessment methods complicates aligning treatment plans with risk of progression. The objective of the study was to evaluate the performance of the 40-GEP test in identifying BWH or AJCC T1 staged tumors.

## Results

**Figure 3. The 40-GEP accurately classifies patients by metastatic risk within AJCC8 and BWH T1 staged subset**



## References


- Karia, et al. JAAD 2012
- Rogers, et al. JAMA Dermatol. 2015
- Alam, et al J Am Acad Dermatol 2018
- Dessinioti et al. J Eur Acad Dermatol Venereol, 2022
- NCCN Guidelines®, v1.2023
- Wysong, et al. J Am Acad Dermatol. 2021
- Ibrahim, et al. Future Oncol. 2021
- Arron, et al. J Drugs Dermatol. 2022
- Farberg, et al. Dermatol Ther. 2022
- Eggermont, et al. J Eur Acad Dermatol, 2023

## Disclosures

› JJS, ALF, and AP are employees and shareholders of Castle Biosciences, Inc.

## Results

**Figure 4. 40-GEP identified ~70% of the T1 tumors that metastasized as biologically high risk**



**Table 1. Distribution of NCCN high risk factors missed by AJCC8/BWH T1 staging that are captured by 40-GEP in real world data**

Risk Factor	% of Patients
Located on Areas H or M*	86.7
Rapidly growing tumor	28.3
Borders poorly defined	17.1
Immuno-suppressed	6.0
Neurological symptoms in tumor region	2.0
Other†	1.7
Tumor at site of prior radiation therapy or chronic inflammatory process	1.2

Test Result	% of Tests	# of Patients
Class 1	80.0%	2,338
Class 2A	18.6%	545
Class 2B	1.4%	41

Clinicopathologic risk factors of patients tested with the 40-GEP that are BWH and AJCC8 T1 (n=2,924) samples. \*mask areas of face, genitalia, hands, feet, cheeks, forehead, scalp, neck and pretibial. † high risk subtypes and infiltrating tumors that are deemed eligible for testing

## Conclusions

- › This study shows that ~70% of the metastasis events that are missed by AJCC8 and BWH T1 staging are captured by 40-GEP Class 2 test results.
- › Real world clinical testing data shows that 20% of patients were staged as T1 and received Class 2 results, identifying potentially 584 patients that could metastasize if only staging criteria were applied and used to determine treatment decisions.
- › These findings demonstrate the ability of the 40-GEP to identify biologically high-risk tumors in patients deemed low risk by traditional staging methods.
- › Incorporating 40-GEP test results in clinical assessments with traditional clinicopathological risk factors can improve the stratification of high-risk cSCC T1 patients and contribute to risk-appropriate surveillance and treatment decisions.